

Charles River Diagnostics through TAGCenter for Transnetyx Users

NOTE: The captured screenshots on this page are of the steps below.

- 1** Place your order in the Charles River ordering platform, Laboratory Testing Management (LTM). Please follow instructions on criver.com/TAGCenter.
- 2** Go to your TAGCenter and sign in via your MyTAG, phone number, or email address
- 3** Scan your Charles River Cell Line Kit as you normally would a Transnetyx well plate.
Note: you have successfully scanned your Charles River Kit when you see the CRL**** barcode appear on the screen.
- 4** Drop your Charles River Cell Line Kit in the TAGCenter slot.
- 5** Confirm your Cell Line Kit was dropped via the confirmation screen.
- 6** Your kit will be shipped by the Central Shippers as normal.

